



Rapid and on-site phytoplankton assessment in water

Sixsenso Technologies S.L

Algae bloom in bathing waters

22 cell/100ml



10 Million cells/100ml



SIXSENSO FOR MONITORING AND PREDICT ALGAE BLOOM OF URTICANT SPECIES LIKE INCREASING OSTREOPSIS IN MEDITERRAENEAN COASTS.



PHYTOPLANKTON SAMPLING SYSTEM: THE WHISH LIST



EASY OF USE

PORTABLE

RAPID

REPRESENTATIVE SAMPLING

EFFICIENT SAMPLING & CONCENTRATION

EASY OF USE

On ship use in remote locations,
automated operation.

PORTABILITY

Small size, Low foot-print and
weight, battery operated.

RAPID

1 liter filtration in 10 minutes.
Powered by 25Watts pump.

REPRESENTATIVENESS

Large Volume (Liters) for low levels
of phytoplankton assessment

EFFICIENCY

Efficient cell filtration and
recovery elution > 95%

SX-CON: SMART SAMPLER CONCENTRATOR

FIRST SIXSENSO'S PRODUCT IN THE MARKET

PRODUCT LAUNCHED IN APRIL 2022



SX-CON is a microorganism concentrator module **which can rapidly filtrate large volume of water samples (in the order of litres) and elute them in a few millilitres reservoir**, providing an increasing concentration factor of some orders of magnitude, depending on the application. **It is designed to be portable (with large battery autonomy) and robust to use in the most challenging outdoors' environmental conditions or simply operating in the laboratory.**

Combine SX-CON with several microbiological assessment techniques and methods in order to decrease time to result, increase their representativeness and enhance their sensitivity, including:

- Lateral Flow Device (LFD) assays,
- Defined Substrate Technology (DST) based assays, Nanoparticles-based assays,
- Antigen-based assays,
- rRNA-based assays,
- culture and
- PCR assays.

SUPPORTED BY



**Agència Catalana
de l'Aigua**

Portable, efficient, robust and easy to use.

"SIXSENSO information for evaluation promotional purposes"





Specifications

Filtrated sample Volume	<i>From 100 mls to 3.000 mls</i>
Eluted sample Volume	<i>From 0,1 ml to 3 mls</i>
Concentration Factor (CF)	<i>From x10 to x1.000</i>
Flow rate	<i>200 ml/min @100 rpm. (Powered by 24W peristaltic pump)</i>
Internal Diameter (ID) tubing	<i>From ID 2.4 to 8.0 mm (ID 8mm by default)</i>
Power	<i>Internal rechargeable battery – 8 hours continuous operation</i>
Dimensions	<i>26,5 x 24,5 x 17 cm</i>
Mass	<i>3 Kg</i>
IP rating lid closed/open	<i>IP68 / IP65</i>

Concentration Kit – Sterilized disposables *

	<i>2x Luer lock 3-port valves</i>
	<i>1x tubing kit</i>
	<i>2x Luer lock 5 ml syringes (with elution buffer)</i>
	<i>1x Celltrap filter</i>



* Sterilized disposables are for only one use to avoid cross-contamination. Multi-use of disposables could be suitable for evaluation purposes.



SX-CYT




Specifications

Laser pump	Diode laser @ 488nm
Detection channels	Green channel @537nm Red channel @ >610nm
Microorganism detection	Bacteria, virus, phytoplankton.
Microorganism sizes	From 0.2 μm to 50 μm .
Operational concentration range	From 1 to 10^4 cell/ml From 1 to 10^4 cell/100ml (with SXS-CON)
Power	Internal rechargeable USB battery – 8 hours continuous operation
Dimensions	36 x 29 x 16 cm
Mass	4 Kg
IP rating lid closed/open	IP68 / IP65



Portable, accurate, robust and easy to use.



- **On-site** microorganism quantification.
- Suitable for **bacteria, virus** and **phytoplankton**.
- **Large sample volume** analysis.
- **Rapid analysis**.
- **Accurate** as gold standard methods.

Phytoplankton: Sampling and Concentration results

PERFORMANCE OF THE CONCENTRATION CONCENTRATING PHYTOPLANKTON FROM REAL BATHING WATERS



1 Liter water sample

Concentration
=
Filtration + Elution



<6 mls concentrated sample providing a Concentration factor in the order of x250.

Optimized concentration protocol has been developed that filters 1L of sample to be further eluted to a volume of ~5.75mL providing a concentration factor in average of **C-Factor: x255 for phytoplankton** with samples up to 10NTU turbidity, with whole process taking **less than 10 minutes**.

“SIXSENSO information for evaluation promotional purposes”


Validation: Phytoplankton concentration & detection









PERFORMANCE OF THE CONCENTRATION AND READER MEASURING PHYTOPLANKTON <50µm IN REAL SEA WATER

Concentration Operation **Range from 10^1 to 10^5 cells/mL** .


- LOD Phytoplankton sub 50 micron:
- **LOD: 53 cells/mL** (95% CI) direct measurement
 - **LOD: 22 cells /100 ml** (95% CI) in combination with the concentrator



Sample of *Tetraselmis* and *Nanocloropsis* phytoplankton shinning under UV light.





Cells of *Tetraselmis* and *Nanocloropsis* phytoplankton passing through the interrogation Field Of View (FOV) of the microfluidic channel of the SX-CYT reader. Samples on the right under UV light showing the eluted concentrated sample shining in red color.



Validation: Environmental Phytoplankton monitoring











Benchmarking against gold standard Flow Cytometry measuring phytoplankton in a Ballast Water Testing ship in the North Sea. Third party validation by MEA-NL



Paper presented at 6th IMarEST Ballast Water Technology Conference at IMO in London (2017).

Validation: Environmental Phytoplankton monitoring

Performance measuring phytoplankton in Ballast Water Vessel discharge Testing in the North Sea.
Third party validation by MEA-NL



"SIXSENSO information for evaluation promotional purposes"

Phytoplankton assessment results presented at the **6th IMarEST Ballast Water Technology Conference at IMO in London (2017)**. Final concentration factor for phytoplankton was x 130, and the **SX-CYT Reader system LoD moved down to 77 cell/L (7,7 cell/100ml)** after the concentration procedure. Achieving a LoD <10 cell/100ml in deep ocean water and meeting with lower alarm level of toxic phytoplankton (DSP < 50 cell/100ml).

References:

J. Pérez, M. Jofre, P.Martinez, A.Paerker, A. Yáñez, V.Catalan, M. Veldhuis, V. Pruneri(2016). "CMOS based image cytometry for detection of phytoplankton in ballast water". [Biomedical Optics Express. 8. 10.1364/BOE.8.001240.](https://doi.org/10.1364/BOE.8.001240)

Pérez, J. M., Jofre, M., Martínez, P., Yáñez, M. A., Catalan, V., & Pruneri, V. (2015). "An image cytometer based on angular spatial frequency processing and its validation for rapid detection and quantification of waterborne microorganisms". *The Analyst*, 140(22), 7734–7741. <https://doi.org/10.1039/c5an01388k>

SIXSENSO PRODUCTS

FOR PREVENTIVE ANALYSIS OF PATHOGENS AND INVASIVE SPECIES
IN MARINE AND FRESH WATERS


- REPRESENTATIVE SAMPLING (LITERS)
- HIGH SENSITIVITY - LOW LIMIT OF DETECTION
- HIGH SPECIFICITY (DNA/RNA LABELLING)
- RAPID ANALYSIS (30 MIN- 3 HOURS)
- IN-SITU / PORTABILITY



PHYTOPLANKTON



BACTERIA



VIRUS



Rapid indicative tests

LOW END PRODUCTS (PASS/FAIL Test):



Filters and tubing
in **SX-BioKit-LFD**



Lateral Flow
Device (LFD) in
the **SX-BioKit-
LFD**



SX-CON

**QUALITATIVE
MICROBIOLOGICAL ASSESSMENT**
< 30 MIN

Quantitative Monitoring System

HIGH END PRODUCTS (Report in cells/100ml)



GENETIC
LABELLING

Reagents and filters in the **SX-BioKit-SYS**



Smart Sampler
Concentrator
Filtration - Elution
SX-CON



Fluorescence
Reader
SX-CYT

**QUANTITATIVE
MICROBIOLOGICAL ASSESSMENT**
< 3 HOURS

SIXSENSO PRODUCTS

AVAILABLE FILTERS FOR SPECIFIC APPLICATION NEEDS

Filtration volumes from 0.1 to 10.000 liters



CT - Filter




CTHV - Filter



Product	CT221	CT401	CTHV 100	CTHV 400	CTHV 1000	CT222	CT402	CT802	CTHV 1002	CTHV 4002	CTHV 10002
Product identifying cap colours	Light blue 	Dark blue 	Yellow 	White 	Dark blue 	Colourless	Red 	White 	Red 	Light blue 	Colourless
Membrane type	PP	PP	PP	PP	PP	PES	PES	PES	PES	PES	PES
Volume (L)	1-5	5-10	1 -100	1 -1,000	1-10,000	0-1	1-5	1-10	1-100	1-1,000	1-10,000
Effective filtration area (cm ²)	20	39	217	815	1,600	11	20	39	108	407	1,080
Hold up volume (mL)	0.4	0.75	3.5	13.3	27	0.1	0.19	0.38	0.89	3.4	8.9
Flow (mL/min)	140	270	1,200	4,482	9,000	242	440	858	2,380	8,950	18,000


IMPACT SCHEME

**BALLAST
WATER**



All is linked. We can **protect food sources, health people, marine biodiversity and climate** by **better monitoring** and controlling invasive alien species (IAS) in Ballast Water operations.



Bacteria, virus and phytoplankton monitoring by Sixsenso systems.



BIODIVERSITY



SEA FOOD



*INVASIVE ALIEN SPECIES (IAS)

**WATER-RELATED INFECTIOUS DISEASES (WRID)

Ballast Water Regulation

REGULATION
ENTERED INTO FORCE
GLOBALLY ON SEPT 2017:
**IMO-BALLAST WATER
D-2 STANDARD**

**TO CONTROL THE INVASIVE SPECIES
BALLASTED TO THE PORT WATER
FROM SHIPS.**



- PHYTOPLANKTON < 10 cell / ml
- E.COLI < 250 CFU / 100ml
- ENTEROCOCCI < 100 CFU/ 100ml
- VIBRIO CHOLERAEE < 1 CFU/100mL

THE TEAM

MANAGEMENT DIRECTORS



Adrian J. Parker | *Co-founder and CEO*

Specialist sales and marketing in water treatment markets.



Pedro Martinez | *Co-founder and CTO*

Photonics, electronics and software development.


TECHNOLOGY ADVISORY BOARD



Prof. Valerio Pruneri |
Co-founder & Technology Advisor
Group Leader @ ICFO

COLLABORATORS

TECHNICAL SUPPORT



Alfredo Ongaro |
Microfluidics

Assay development and microfluidics design



David Kernan | *RD*

Protocol development for sample filtration and concentration



PUBLIC SUPPORT



Agència Catalana de l'Aigua



ACCIÓ



PRIVATE ALLIANCES





THANKS